





TELEDYNE CETAC TECHNOLOGIES
Everywhereyoulook™

Liquid Sample Introduction for ICP-AES

Measurement of Elemental Impurities in a Higher Daily Dose Drug Product by USP <232>/<233> using Ultrasonic Nebulization with ICP-AES Detection



USP <232>/<233>

On January 1, 2018 the U.S. Pharmacopeia (USP) enacted new criteria for element impurities in finished drug products. These criteria, detailed in USP <232>/<233>, recommend analysis of drug products for element impurities by either Inductively Coupled Plasma Atomic Emission Spectroscopy (ICP-AES) or Inductively Coupled Plasma Mass Spectrometry (ICP-MS).

Laboratories must measure impurities based on a J value for each drug product; drug products may be oral, parenteral (ex. intravenous, injection), or inhalation types. The J value of the drug is calculated based on an established permissible daily exposure (PDE), maximum daily dose (MDD) of the drug, and the dilution factor used in the sample preparation method. As a result of the J value calculation, drugs with a larger MDD require lower element impurity detection limits.

Liquid Sample Introduction for ICP-AES

Ultrasonic Nebulization

An ultrasonic nebulizer is an accessory for ICP-AES that enables higher sample transport efficiency (versus a standard pneumatic nebulizer) to the ICP-AES plasma. This benefit can be helpful for detection of more difficult elements such as As, Cd, Pb, Hg, and Tl.

This work describes the use of ultrasonic nebulization for ICP-AES detection of element impurities in aspirin, a drug with a higher daily dose than a low dose drug (ex. 1 tablet per day) such as an allergy or sleep aid medicine.

Liquid Sample Introduction

Ultrasonic Nebulization

In place of a regulated gas flow for generation of a liquid sample aerosol (pneumatic nebulization), liquid sample is pumped across a quartz plate with an underlying oscillating (piezoelectric) crystal.

The oscillations of the crystal will break up the liquid flow and cause formation of a sample aerosol. Ultrasonic nebulization is typically up to 10x more efficient (versus a conventional pneumatic nebulizer) for conversion of liquid sample into a useable aerosol.

U5000AT⁺ Ultrasonic Nebulizer

Dimensions and Weight



Width: 35.6 cm (14.0 in)

Depth: 34.9 cm (13.7 in)

Height: 25.4 cm (10.0 in)


Weight: 12.3 kg (27 lbs)

 **TELEDYNE**
CETAC TECHNOLOGIES
Ultrasonic Nebulizer



Teledyne CETAC U5000AT+ Ultrasonic Nebulizer

Front View



**Temperature
Controllers**

Teledyne CETAC U5000AT+ Ultrasonic Nebulizer

Back View

Sample
Out


Argon
In

Power
Input

Drain
Pump




Teledyne CETAC U5000AT+ Schematic



**Piezoelectric Crystal
(Transducer)**


U5000AT+ Transducer Assembly



**Piezoelectric
Crystal
(disk shape)**


Quartz Plate

U5000AT⁺ Ultrasonic Nebulizer Desolvation System



U5000AT⁺ Ultrasonic Nebulizer

Example Interface Kit



U5000AT⁺ Ultrasonic Nebulizer

Installation Steps

1. Remove standard nebulizer and spray chamber from the ICP-AES.
2. Connect nebulizer gas inlet line from ICP-AES to the USN.
3. Connect sample out line from USN to ICP torch.
4. Connect sample line from autosampler to USN inlet line.
5. Connect power cord to USN and turn on. Wait 10 minutes for heater and cooler to stabilize. Start ICP, introduce tune solution, press one button (Operate) to begin nebulizing samples.

Overall setup takes about 15 minutes.


No computer control, no software installation needed.

ICP-AES w. Ultrasonic Nebulization



Agilent 5100 ICP-AES & U5000AT+ USN

ICP-AES w. Ultrasonic Nebulization




ICP-AES w. Ultrasonic Nebulization




PerkinElmer Avio 500 ICP-AES, U5000AT+ USN & ASX-280 Autosampler

Connections Between Host ICP-AES and U5000AT+ Ultrasonic Nebulizer



Example Torch Adapter Connection



12/5 Glass
Socket
Adapter

Instrumentation

ICP-AES: PerkinElmer Avio 500

Ultrasonic Nebulizer (USN): Teledyne CETAC U5000AT⁺


Microwave Digestion System: CEM Mars 6

U5000AT⁺ Ultrasonic Nebulizer

Aerosol Generation



Ultrasonic Nebulizer Aerosol Close-Up View



USP <232> Elements Class Grouping

- **Class 1 and 2A elements must be measured; Class 1 elements are most toxic (Cd, Pb, As, Hg).**
- **Class 2B elements must be measured if added to the drug product generation process.**
- **Class 3 elements are of lower toxicity via oral administration, but require measurement if drug product given by inhalation or parenteral routes.**

Maximum Oral Daily Exposures for Elements Defined in USP <232> - I

Element	Class	Max. Oral Daily Exposure ($\mu\text{g}/\text{day}$)
Cd	1	5
Pb	1	5
As	1	15
Hg	1	30
Co	2A	50
V	2A	100
Ni	2A	200
Tl	2B	8
Au	2B	100
Pd	2B	100
Ir	2B	100
Rh	2B	100

Maximum Oral Daily Exposures for Elements Defined in USP <232> - II

Element	Class	Max. Oral Daily Exposure ($\mu\text{g}/\text{day}$)
Ru	2B	100
Se	2B	150
Ag	2B	150
Pt	2B	100
Li	3	550
Sb	3	1200
Ba	3	1400
Mo	3	3000
Cu	3	3000
Sn	3	6000
Cr	3	11000

Sample Preparation

A 0.5 g aspirin sample with multielement spike was added to each digestion vessel followed by 5 mL of reagent grade HNO₃ and 1 mL of reagent grade HCl. The vessels were left uncapped for 10 minutes in a fume hood to allow any initial gases to vent prior to sealing the vessels. Spiked and unspiked samples were digested using a closed vessel microwave digestion program as specified by USP <233>.

Stage	Power(W)	Ramp	Hold	Temp (°C)
1	1050	15 min	15 min	200

Sample Dilution and J Values

The final dilution for each sample after digestion was 100x with deionized water. With that dilution factor and a maximum daily dose of 4.32 g (one aspirin tablet is 0.36g, maximum dose is 12 tablets per day), the calculated J values (rounded down) are listed in the next two slides in mg/L. Following USP <233> protocol, a reagent blank, 0.5 J standard, and 1.5 J standard were used for calibration.

Analytical Concentrations at Different J Values - I

Element	J-value (mg/L)	0.5 J	1.5 J
Cd	0.01	0.005	0.015
Pb	0.01	0.005	0.015
As*	0.03	0.015	0.06
Hg*	0.06	0.03	0.09
Co	0.10	0.05	0.15
V	0.20	0.10	0.30
Ni	0.40	0.20	0.60
Tl	0.016	0.008	0.024
Au	0.20	0.10	0.30
Pd	0.20	0.10	0.30
Ir	0.20	0.10	0.30
Rh	0.20	0.10	0.30

Analytical Concentrations at Different J Values - II

Element	J-value (mg/L)	0.5 J	1.5 J
Ru	0.20	0.10	0.30
Se	0.30	0.15	0.45
Ag	0.30	0.15	0.45
Pt	0.20	0.10	0.30
Li	1.1	0.55	1.65
Sb	2.4	1.2	3.6
Ba	2.8	1.4	4.2
Mo	6	3	9
Cu	6	3	9
Sn	12	6	18
Cr	22	11	33


Calibration & Reagent Addition

The ICP-AES was calibrated using standards that were matrix matched to the acid concentrations of the digested samples. An internal standard solution of 50 µg/L Y was added to all samples and standards using a mixing tee.

Samples for the Class 1 elements (includes Hg), 2A elements, and thallium were digested and analyzed separately from the Class 2B and Class 3 elements. For these elements, the reagent **L-cysteine, a thiol-containing amino acid**, was added to digested samples to assist element transport of Hg through the ultrasonic nebulizer.

L-cysteine was added to the 50 µg/L Y internal standard solution so the final concentration after sample mixing was 3 mg/ml. As nitric acid rapidly oxidizes L-cysteine, the internal standard solution was prepared in 0.07M HCl.

L-Cysteine (C₃H₇NO₂S) Structure



Thiol functional group 

Aspirin Tablets: Notes

- Aspirin tablets did not contain SiO_2 or TiO_2 , so hydrofluoric acid (HF) not needed for digestion
- If HF was required for digestion, then a neutralization of residual HF would be needed (ex. boric acid) before introduction to any glass components

Emission Line Selection

- All elements required by USP <232> and other elements potentially in the aspirin sample (ex. Ca, Na, Si, Mg) were run individually to check for spectral interferences.
- When possible emission lines free of interferences were selected; otherwise multicomponent spectral fitting (MSF) or interelement correction factors (IECs) were used.
- An interference check solution was measured during the sample run to verify that spectral correction techniques were working properly.

Operating Conditions

ICP-AES with Standard Nebulizer and U5000AT+ USN

Parameter	Std Nebulizer	U5000AT+ USN
ICP Power	1500 W	1500 W
Plasma Gas	8.0 L/min	8.0 L/min
Auxiliary Gas	0.2 L/min	0.2 L/min
Nebulizer Gas	0.70 L/min	0.62 L/min
Torch Injector	2 mm	2 mm
Uptake Rate	1.5 mL/min	1.0 mL/min
Cassette Position	-3.0	-5.0
Resolution	Normal	Normal
Nebulizer Type	Meinhard K	Piezoelectric
Spray Chamber	Baffled cyclonic	Conical
Heater Temp	N/A	120°C
Cooler Temp	N/A	5°C
Integration Time	2 s min, 10 s max	2 s min, 10 s max
Peak Area	3 pts/peak	3 pts/peak
Replicates	3	3


Elements, Wavelengths, Plasma View - I

Element	Wavelength (nm)	Plasma View
Cd	228.802	Axial
Pb	220.353	Axial
As	188.979	Axial
Hg	253.652	Axial
Co	228.616	Axial
V	292.402	Axial
Ni	231.604	Axial
Tl	190.801	Axial
Au	242.795	Axial
Pd	340.458	Radial
Ir	208.882	Axial
Rh	343.489	Radial


Elements, Wavelengths, Plasma View - II

Element	Wavelength (nm)	Plasma View
Ru	349.894	Radial
Se	196.026	Axial
Ag	338.289	Axial
Pt	265.945	Axial
Li	670.784	Radial
Sb	231.146	Radial
Ba	493.408	Radial
Mo	202.030	Axial
Cu	324.752	Radial
Sn	189.927	Axial
Cr	267.716	Radial


Pb (220.353 nm) Spectra at 0.01 mg/L Pb (Class 1)




As (188.979 nm) Spectra at 0.03 mg/L As (Class 1)



Hg (253.652 nm) Spectra at 0.06 mg/L Hg (Class 1)



TI (190.801 nm) Spectra at 0.016 mg/L TI (Class 2B)




--- USN
--- Std Neb

Analytical Criteria Defined in USP <233> for Quantitative Procedures


Criteria	Description
Accuracy	Spike recoveries at 0.5J, J, and 1.5J must be between 70-150%
Repeatability	The %RSDs of measurements of six independent samples spiked at J must be less than 20%
Ruggedness	Six solutions must be analyzed on different days, with different instruments, or with different analysts. The %RSDs over the 12 measurements must be less than 25%
System Suitability	The difference in the results of the high calibration standard (1.5J) measured at the beginning and end of a batch must be < 20%

Accuracy: 0.5J and 1.5J Spike Recoveries in Aspirin


Black (top) and Gray (bottom) Lines Show USP <233> Limits




Repeatability: %RSDs of Six Independent Samples Spiked at 1J



Ruggedness: %RSDs of 6 Independent Samples Analyzed Over 2 Days (12 Measurements)



System Suitability: Analysis of 1.5J at the Beginning and End of Analytical Run



Difference in the results must be < 20%.

Summary

All criteria required by USP <233> were met using the U5000AT⁺ Ultrasonic Nebulizer with ICP-AES detection of trace elements in aspirin. The U5000AT⁺ enables enhanced analyte sensitivity and lower background emission, improving ICP-AES measurement of more difficult elements in this example of a higher daily dose drug product.

U5000AT+ USN: Benefits & Advantages

- **Enhanced detection of difficult Class 1 elements such as Pb, As, Hg and Class 2B element Tl**
- **Lower background emission due to less injected water vapor**
- **Fast and easy setup of the U5000AT+**
- **Convenient placement on laboratory cart or benchtop**
- **One button operation**
- **No computer control required**
- **Simple connection to autosampler (ex. ASX-280/560) for automated sample introduction**

Ultrasonic Nebulizers for ICP-AES

Top Level Part Numbers

- **U5000AT⁺ Ultrasonic Nebulizer**
 - **U51-99-0001A+ (115 V)**
 - **U51-99-0001B+ (220 V)**

U5000AT+ Interface Kits for ICP-AES

U5000AT+ Compatibility

- **All ICP-AES Models**
 - **Examples with Interface Kit SP Numbers ()**
 - **Agilent (SP5155Y)**
 - **Analytik Jena (SP5155Y)**
 - **PerkinElmer (SP5155L)**
 - **Shimadzu (SP5155Z)**
 - **Spectro (SP5155P)**
 - **Teledyne Leeman (SP5155I)**
 - **Thermo Fisher (SP5155CC)**

One interface kit provided at no charge with each U5000AT+

Acknowledgement

- **Special thanks to Mary Jo Menke-Wright for her work in sample preparation, method development, data collection, and outline text for this important application.**

Where to go for more information

- Service support
 - cetacservice@teledyne.com
- Sales support
 - cetacsales@teledyne.com
- Technical support
 - Fred.Smith@Teledyne.com
 - Paula.Doeschot@Teledyne.com
- Web information
 - <http://www.teledynecetac.com/products/nebulizers/u5000at+>
 - <http://www.teledynecetac.com/resourceSite/Application%20Notes/AP-U5000-003.pdf>